

Health Alerts

The U.S. Consumer Product Safety Commission (CPSC) has announced the recall of the following products. Consumers should stop using recalled products unless otherwise instructed. Consumers can submit reports of harm to CPSC's new searchable online product safety database at www.SaferProducts.gov. A searchable food and medical product recall database is available at www.fda.gov/Safety/Recalls/default.htm.

Breast pump power adapters

Hazard: The casing on certain AC/DC power adapters that are used with the Playtex Nurser Deluxe Double Electric Breast Pump may become loose and separate, causing an electric shock.

Description: The product can be identified by product serial number (P12324-XXXX through P13205-XXXX) or by the adapter production code (1241 through 1324).

Remedy: Contact Playtex for a replacement.

Contact: Call 888-207-1492 or visit www.playtexproducts.com/bppoweradapter.aspx.

